

## CMPT Mk 2 Positioning Options

1



The Trainer can be used 2 ways:

- standing
- supine (using the supplied base)

## Clinical Male Pelvic Trainer (CMPT) Mk 2 - Standard



Designed and manufactured by Limbs & Things Limited,  
Sussex Street, St. Philips, Bristol, BS2 0RA, UK.  
Telephone: +44 (0)117 311 0500 Fax: +44 (0)117 311 0501  
sales@limbsandthings.com www.limbsandthings.com

Part No: **160-000**

Issue 2, September 2014  
© 2014 Limbs & Things

For more skills training products visit

**[www.limbsandthings.com](http://www.limbsandthings.com)**

# Clinical Male Pelvic Trainer Mk 2

## Standard Part No: 60950

The Clinical Male Pelvic Trainer Mk 2 presents key anatomical features, both externally and internally, for teaching 'hands-on' male pelvic examination and diagnosis. The uncluttered design of the model enables trainees to focus on essential anatomy and procedure. The Mk 2 model features new and improved pathologies which allow for a realistic and anatomically accurate training experience. This product has been designed in close collaboration with the staff at the Clinical Skills Resource Centre, University of Liverpool; Imperial College, London; Guys Hospital, London; Southmead Hospital, Bristol; University of Southampton, UK.

### Skills

- Learning examination procedure
- Testicular examination
- Dry catheterisation (Foley catheter size 16 and retrograde catheters)
- Examination of normal anatomy
- Evaluate testicular abnormalities
- Evaluate abdominal / pelvic pain

### Features

- The genitalia includes the penis (both circumcised and non-circumcised), scrotum, testicles (includes vas deferens and epididymis)
- Anatomy of penis features glans and shaft
- Anatomy of groin with clear anatomical landmarks features ASIS (anterior superior iliac spine) and pubic symphysis with tubercles
- Anatomy of scrotum area: scrotum; testicles with epididymis and vas deferens
- Can be used in two positions (standing and supine)
- Soft tissue inserts are removable and replaceable: genitals (penis, perineum, & testicles) and abdominal wall
- Skin surface is washable using soap and water

## Package supplied Part No

### Standard

- 1 **CMPT Module 1 Normal** 60952
- 1 **CMPT Abdominal Insert** 60959
- 1 **CMPT Pubic Bone** 60960
- 2 **CMPT Foreskins (Pack of 2)** 60961
- 1 lower torso
- 1 base (for supine use)
- 3 sachets of waterbased lubricant
- 1 carry case

## Components



**CMPT Abdominal Insert**

Part No: **60959**



**CMPT Module 1: Normal**

Part No: **60952**



**Lower torso**



**Base**



**CMPT Foreskins (Pack of 2)**

Part No: **60961**



**CMPT Pubic Bone**

Part No: **60960**


## Initial Setup

During shipping the genital module is not fully attached, hence before first use, please ensure you hook the top skin flap onto the pubic bone locating pin.

1




2



## Changing the CMPT Modules

### 1 Removing



Remove the Abdominal Insert.

Hold the Module and Pubic Bone assembly at the front and gently remove from the lower torso.

### 2 Refitting



When refitting, make sure that the Pubic Bone is properly seated in the recess in the lower torso.

When refitting the Abdominal Insert fit the skin over the the back edge of the torso first and then lower the Insert down into position.


### 3



The Trainer is ready to use.



## Removing/Attaching CMPT Module to Pubic Bone

### 1



Lift the skin flap on the back of the CMPT Module away from the locating pin on the Pubic Bone.


### 2



Gently pull the Module and Pubic Bone apart.

To reattach a Module to the Pubic Bone, push the pubic arch down onto the inverted 'V' of the Module.

### 3





Ensure that the Pubic Bone is secured behind the lip of the V on the Module.

Finally pull the skin flap back onto the locating pin on the Pubic Bone.



# Attaching the Foreskin

1



Using the two main fingers of each hand, stretch the end of the Foreskin so that it forms a square approximately 5cm across.

2



Keeping the Foreskin stretched, line up the hole in the Foreskin with the end of the penis. Now pull the Foreskin just past the glans.

3





The narrowest part of the Foreskin should be at the tip of the penis.

Add some lubricant to the end of the penis.


**Note: Do not rub the lubricant around the end of the penis as this will make handling the Foreskin difficult in the next stage.**

4



Pull the open end of the foreskin back down over the length of the penis.

5



The end of the Foreskin should finish at the base of the penis. Wipe any excess lubricant from the base of the penis.

Massage the foreskin to distribute the lubricant thoroughly around it.

6



The foreskin can now be smoothly and realistically retracted.

**Note: Do not leave the Foreskin on the model as the waterbased lubricant will dry out over time and stop it functioning properly.**